


InnoVESP™


Disposable Endocavity Needle Guides




Redefine **Ovum Pick-up**

InnoFine

Innovative design




Pure metal track




Visualized

- Ultra wide view for needle insertion




Efficient

- Faster operations
- Easier detach after procedures



Safe

- No microplastics into human body
- Sing-use design reduces risk of cross-contamination



Precise

- To the exact region of interest

InnoVESP™

28 models for Total Solution

REF#	Compatible transducers	Needle size/ Channel length
DE-001	【Mindray】 6CV1, V10-4, V11-3E, V11-3HU 【Samsung】 EC4-9IS 【Wisonic】 EV10-4 【VINNO】 G4-9E 【Landwind】 EC4-9	16-18G/15cm
DE-002	【GE】 E8C, E8CS, IC5-9-D 【Samsung】 VR5-9	16-18G/15cm
DE-003	【Philips】 C8-4v, C10-3v, C9-4v(ClearVue)	16-18G/15cm
DE-004	【GE】 IC9-RS	16-18G/15cm
DE-005	【Hitachi-Aloka】 UST-984-5, UST-981-5, C41B	16-18G/21.5cm
DE-006	【Esaote】 EC123, EC1123(3.8°)	16-18G/15cm
DE-007	【Esaote】 EC1123(0°)	16-18G/15cm
DE-008	【Canon/Toshiba】 661VT 【Siemens】 EC9-4(ACUSON S1000, S2000, S3000 and Sonoline Antares)	16-18G/15cm
DE-009	【Hitachi-Aloka】 EUP-V53W, EUP-CC531, C41V, CC41R	16-18G/15cm
DE-010	【Hitachi-Aloka】 EUP-V73W, C41V1	16-18G/16.5cm
DE-011	【Philips】 C9-4v(Affiniti, SPARQ)	16-18G/15cm
DE-012	【Samsung】 EVN4-9, ER4-9 【ALPINION】 EN3-10 【Ultrasonix】 EC9-5/10	16-18G/15cm
DE-013	【Canon/Toshiba】 781VT, 781VTE	16-18G/15cm
DE-014	【Konica-Minolta】 EC9-3 【Siemens】 MC9-4, 9MC3, 9EC4	16-18G/15.5cm
DE-015	【GE】 RIC5-9-D, RIC5-9A-RS 【Samsung】 V5-9 【Mindray】 DE10-3U	16-18G/16cm
DE-016	【SonoScape】 6V7	16-18G/21.5cm
DE-017	【Philips】 BP10-5ec 【Samsung】 E3-12A, NEV4-9ES 【ALPINION】 E3-10H	16-18G/15.5cm
DE-018	【Hitachi-Aloka】 UST-9118, UST-9124	16-18G/16cm
DE-019	【SonoScape】 6V3	16-18G/16cm
DE-020	【SuperSonic】 SE12-3	16-18G/15.5cm
DE-021	【Samsung】 EA2-11B 【Siemens】 EC9-4(ACUSON X150 and X300) 【Terason】 8EC4A	16-18G/15cm
DE-022	【Philips】 C9-5ec	16-18G/11cm
DE-023	【Philips】 C10-4ec	16-18G/11cm
DE-024	【WELLD】 8EC4A	16-18G/15cm
DE-025	【SIUI】 V6L11C	16-18G/15.5cm
DE-026	【EDAN】 E10-3B	16-18G/21.5cm
DE-027	【ALPINION】 EV3-10 【Neusoft】 EV3-10X	16-18G/16cm
DE-028	【Siemens】 EC10-5w	16-18G/16.5cm



Leading Provider of Biopsy Needle Guides

InnoFine

JingFang Precision Medical Device(Shenzhen)Co.,Ltd.

3/F,Bldg C,JinWeiYuan Industry Park,Qinglan 2nd Road,
Pingshan,Shenzhen 518118,China

T : 0086(0)0755 8945 8526

E : sales@innofine.cn

w : www.innofine.cn