

For different sizes of interventional instruments, 10-23G/6-10F

Supporting accurate and angle adjustable biopsy

Needle Guides (Disposable)

Bracket (Reusable)

Elastic buckle for safe attaching to the transducer

Simple pressing quick release

Disposable Needle Guides for In-plane ultrasound guided interventions

- Liver, kidney biopsy
- Tumor ablation
- Aspiration and drainage of cyst, abscess, effusions
- Thyroid , breast and lymph nodes biopsy
- Anesthesia

Disposable Needle Guides for general purpose biopsy

(JSP is Reusable Bracket , DGC is Disposable Needle Guides)

Bracket REF#	Compatible Transducers	Disposable Needle Guides
JSP-047	【GE Healthcare】 4C,4C-A,4C-D,4C-RS,4C-SC 【Samsung】 CS1-4	<p>DGC-01 for 10-23G/6-10F</p> <p>DGC-0101 for 10G/10F, 12G/8F,14G/6F,16G,18G, 20G,22G</p> <p>DGC-0102 for 11G/9F, 13G/7F,15G,17G,19G,21G, 23G</p> <p>Note: Each JSP reusable bracket with 5 packs of DGC-0101 disposable needle guide</p>
JSP-001	【GE Healthcare】 C1-5-D,C1-5-RS	
JSP-016	【GE Healthcare】 C1-6-D,C1-6VN-D	
JSP-028	【GE Healthcare】 11L,12L,11L-D,12L-RS,12L-SC,L6-12-RS	
JSP-051	【GE Healthcare】 ML6-15-D 【Samsung】 LM4-15B 【Mindray】 LM15-6U	
JSP-050	【GE Healthcare】 C2-7-D,C2-7VN-D	
JSP-018	【Philips】 C5-1(CX50,EPIQ5,EPIQ7)	
JSP-033	【Philips】 L12-4(Affiniti50,Affiniti70 and ClearVue),L12-3(Affiniti50,Affiniti70,EPIQ5, EPIQ7 and CX50)	
JSP-052	【Hitachi-Aloka】 UST-9123,UST-9130	
JSP-053	【Hitachi-Aloka】 C251,C35,C253	
JSP-054	【Hitachi-Aloka】 UST-9133,C22K	
JSP-031	【Canon】 PVT-375AT,PVT-375AX,PVT-375BT,PVT-375ST,PVU-375BT,PVT-375SC	
JSP-055	【Canon】 PVT-382BT,PVU-382BT	
JSP-022	【Mindray】 C5-1s,C5-1E,C5-1U,SC6-1E,SC6-1U	
JSP-007	【Mindray】 L12-4s,L12-3E	
JSP-035	【Samsung】 CA1-7A,CA1-7AD	
JSP-041	【Samsung】 C2-8A,CA2-8A,CA2-8AD,CA2-9A,CA2-9AD,CF2-8	
JSP-029	【Siemens】 4C1,6C1HD 【Samsung】 C1-4EC	
JSP-024	【Siemens】 4V1,PH4-1	
JSP-002	【Esaote】 CA541	
JSP-043	【SONOSCAPE】 3C-A	
JSP-026	【VINNO】 X6-16L	
JSP-039	【VINNO】 X4-12L	
JSP-038	【CHISON】 D7L40L	
JSP-023	【Terason】 5C2-A	
JSP-045	【ALPINION】 SC1-4H,C1-6,SC1-6,SC1-6H 【eZono】 35CL60,CL1-6 【Neusoft】 SC1-6,SC1-6H	
JSP-044	【ALPINION】 L3-12,L3-12T,L3-12i 【eZono】 8LA,L3-12 【Neusoft】 L3-12,L3-12T,L3-12i	

Please contact InnoFine for more models.



JingFang Precision Medical Device(Shenzhen) Co., Ltd.

Address: 3/F, Bldg C, Jinwei Yuan, Crossing at Qingsong West Road and Qinglan 2nd Road, Julongshan, Grand Industrial Zone, Longtian Street, Pingshan District, Shenzhen 518118, China

Tel: +86 (0)755-27597952

Email: sales@innofine.cn

Website: www.innofine.cn